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8. Referências

1. Hilgers FJ, Ackerstaff AH, Aaronson NK, Schouwenburg PF, Van ZN. Physical and psychosocial consequences of total laryngectomy. *Clin Otolaryngol Allied Sci.* 1990;15(5):421-5.
2. McFadden ER, Jr., Denison DM, Waller JF, Assoufi B, Peacock A, Sopwith T. Direct recordings of the temperatures in the tracheobronchial tree in normal man. *J Clin Invest.* 1982;69(3):700-5.
3. Ingelstedt S. Studies on the conditioning of air in the respiratory tract. *Acta Otolaryngol Suppl.* 1956;131:1-80.
4. Keck T, Leiacker R, Heinrich A, Kuhnemann S, Rettinger G. Humidity and temperature profile in the nasal cavity. *Rhinology.* 2000;38(4):167-71.
5. Scheenstra RJ, Muller SH, Vincent A, Hilgers FJ. Heat and moisture exchange capacity of the upper respiratory tract and the effect of tracheotomy breathing on endotracheal climate. *Head Neck.* 2011;33(1):117-24.
6. Zuur JK, Muller SH, de Jongh FH, Van ZN, Hilgers FJ. The physiological rationale of heat and moisture exchangers in post-laryngectomy pulmonary rehabilitation: a review. *Eur Arch Otorhinolaryngol.* 2006;263(1):1-8.
7. Lorenz KJ. *Pulmonary Rehabilitation after Laryngectomy.* 1st edition - Bremen: UNI-MED Verlag AG; 2013 (UNI-MED SCIENCE). 86 p.
8. Williams R, Rankin N, Smith T, Galler D, Seakins P. Relationship between the humidity and temperature of inspired gas and the function of the airway mucosa. *Crit Care Med.* 1996;24(11):1920-9.
9. Mercke U. The influence of varying air humidity on mucociliary activity. *Acta Otolaryngol.* 1975;79(1-2):133-9.
10. Mercke U, Toremalm NG. Air humidity and mucociliary activity. *Ann Otol Rhinol Laryngol.* 1976;85(1 Pt 1):32-7.
11. Schwab JA, Zenkel M. Filtration of particulates in the human nose. *Laryngoscope.* 1998;108(1 Pt 1):120-4.
12. Hatch TF. Distribution and deposition of inhaled particles in respiratory tract. *Bacteriol Rev.* 1961;25:237-40.
13. Tellier R. Aerosol transmission of influenza A virus: a review of new studies. *J R Soc Interface.* 2009;6 Suppl 6:S783-90.
14. Kelly FJ, Fussell JC. Air pollution and airway disease. *Clin Exp Allergy.* 2011;41(8):1059-71.
15. Particulate matter: Basic Information: US EPA (United States Environmental Protection Agency); 2014 [atualizado 2014. Disponível em: <http://www.epa.gov/pm/basic.html>.
16. McGarry P, Morawska L, He C, Jayaratne R, Falk M, Tran Q, et al. Exposure to particles from laser printers operating within office workplaces. *Environ Sci Technol.* 2011;45(15):6444-52.

17. Invernizzi G, Ruprecht A, Mazza R, Rossetti E, Sasco A, Nardini S, et al. Particulate matter from tobacco versus diesel car exhaust: an educational perspective. *Tob Control*. 2004;13(3):219-21.
18. Air quality and health - Fact sheet No 313 2014 [updated 2014. Disponível em: <http://www.who.int/mediacentre/factsheets/fs313/en/>.
19. The World Health Report 2002 - Reducing Risks, Promoting Healthy Life. 2002 [updated 2002. Disponível em: <http://whqlibdoc.who.int/publications/2002/9241562072.pdf>.
20. WHO Air quality guidelines for particulate matter, ozone, nitrogen dioxide and sulfur dioxide, Global update 2005, Summary of risk assessment. 2005 [updated 2005. Disponível em: http://whqlibdoc.who.int/hq/2006/WHO_SDE_PHE_OEH_06.02_eng.pdf.
21. Anderson HR, Spix C, Medina S, Schouten JP, Castellsague J, Rossi G, et al. Air pollution and daily admissions for chronic obstructive pulmonary disease in 6 European cities: results from the APHEA project. *Eur Respir J*. 1997;10(5):1064-71.
22. Brunekreeft B, Holgate ST. Air pollution and health. *Lancet*. 2002;360(9341):1233-42.
23. Jongeneel WP, Straatsen BAM, van Kempen EEMM, Fischer PH. Gezondheidseffecten van wegverkeer: een quickscan 2008 [updated 2008. Disponível em: <http://www.ggd Kennisnet.nl/?file=5785&m=1319011665&action=file.download>.
24. Villeneuve PJ, Chen L, Rowe BH, Coates F. Outdoor air pollution and emergency department visits for asthma among children and adults: a case-crossover study in northern Alberta, Canada. *Environ Health*. 2007;6:40.
25. Lepeule J, Laden F, Dockery D, Schwartz J. Chronic exposure to fine particles and mortality: an extended follow-up of the Harvard Six Cities study from 1974 to 2009. *Environ Health Perspect*. 2012;120(7):965-70.
26. Anderson JO, Thundiyil JG, Stolbach A. Clearing the air: a review of the effects of particulate matter air pollution on human health. *J Med Toxicol*. 2012;8(2):166-75.
27. Heyder J, Rudolf G. Deposition of aerosol particles in the human nose. *Inhaled Part*. 1975;4 Pt 1:107-26.
28. Londahl J, Moller W, Pagels JH, Kreyling WG, Swietlicki E, Schmid O. Measurement techniques for respiratory tract deposition of airborne nanoparticles: a critical review. *J Aerosol Med Pulm Drug Deliv*. 2014;27(4):229-54.
29. Keck T, Leiacker R, Klotz M, Lindemann J, Riechelmann H, Rettinger G. Detection of particles within the nasal airways during respiration. *Eur Arch Otorhinolaryngol*. 2000;257(9):493-7.
30. Jay S, Ruddy J, Cullen RJ. Laryngectomy: the patient's view. *J Laryngol Otol*. 1991;105(11):934-8.
31. van den Boer C, van Harten MC, Hilgers FJ, van den Brekel MW, Retel VP. Incidence of severe tracheobronchitis and pneumonia in laryngectomized patients: a retrospective clinical study and a European-wide survey among head and neck surgeons. *Eur Arch Otorhinolaryngol*. 2014;271(12):3297-303.

32. Retel VP, van den Boer C, Steuten LM, Okla S, Hilgers FJ, van den Brekel MW. Cost-effectiveness of heat and moisture exchangers compared to usual care for pulmonary rehabilitation after total laryngectomy in Poland. *Eur Arch Otorhinolaryngol*. 2015;272(9):2381-8.
33. McRae D, Young P, Hamilton J, Jones A. Raising airway resistance in laryngectomees increases tissue oxygen saturation. *Clin Otolaryngol Allied Sci*. 1996;21(4):366-8.
34. Ferris BG, Jr., Mead J, Opie LH. Partitioning of respiratory flow resistance in man. *J Appl Physiol*. 1964;19:653-8.
35. Recinto C, Efthymeou T, Boffelli PT, Navalta JW. Effects of Nasal or Oral Breathing on Anaerobic Power Output and Metabolic Responses. *Int J Exerc Sci*. 2017;10(4):506-14.
36. Togawa K, Konno A, Hoshino T. A physiologic study on respiratory handicap of the laryngectomized. *Arch Otorhinolaryngol*. 1980;229(1):69-79.
37. Scheenstra RJ, Muller SH, Vincent A, Sinaasappel M, Hilgers FJ. Influence of breathing resistance of heat and moisture exchangers on tracheal climate and breathing pattern in laryngectomized individuals. *Head Neck*. 2010;32(8):1069-78.
38. McRae RD, Jones AS, Young P, Hamilton J. Resistance, humidity and temperature of the tracheal airway. *Clin Otolaryngol Allied Sci*. 1995;20(4):355-6.
39. Swift AC, Campbell IT, McKown TM. Oronasal obstruction, lung volumes, and arterial oxygenation. *Lancet*. 1988;1(8577):73-5.
40. Harris S, Jonson B. Lung function before and after laryngectomy. *Acta Otolaryngol*. 1974;78(3-4):287-94.
41. Natvig K. Influence of different climates on the peak expiratory flow in laryngectomees. *J Laryngol Otol*. 1984;98(1):53-8.
42. Todisco T, Maurizi M, Paludetti G, Dottorini M, Merante F. Laryngeal cancer: long-term follow-up of respiratory functions after laryngectomy. *Respiration*. 1984;45(3):303-15.
43. Torjussen W. Airway obstructions in laryngectomized patients. A spirometric investigation. *Acta Otolaryngol*. 1968;66(1):161-70.
44. Usui N. Ventilatory function in laryngectomized patients. *Auris Nasus Larynx*. 1979;6(2):87-96.
45. Foreman A, De Santis RJ, Sultanov F, Enepekides DJ, Higgins KM. Heat and moisture exchanger use reduces in-hospital complications following total laryngectomy: a case-control study. *J Otolaryngol Head Neck Surg*. 2016;45(1):40.
46. Griffith TE, Friedberg SA. Histologic Changes in the Trachea Following Laryngectomy. *Ann Otol Rhinol Laryngol*. 1964;73:883-92.
47. Roessler F, Grossenbacher R, Walt H. [Condition of the tracheobronchial mucous membrane in patients with long-term tracheostomy. A scanning electron microscopy study]. *Laryngol Rhinol Otol (Stuttg)*. 1988;67(2):66-71.
48. Zuur JK, Muller SH, Sinaasappel M, Hart GA, Van ZN, Hilgers FJ. Influence of heat and moisture exchanger respiratory load on transcutaneous oxygenation in laryngectomized individuals: a randomized crossover study. *Head Neck*. 2007;29(12):1102-10.

49. Zuur JK, Muller SH, Vincent A, Sinaasappel M, de Jongh FH, Hilgers FJ. The influence of a heat and moisture exchanger on tracheal climate in a cold environment. *Med Eng Phys.* 2009;31(7):852-7.
50. Ackerstaff AH, Hilgers FJ, Balm AJ, Van Zandwijk N. Long-term pulmonary function after total laryngectomy. *Clin Otolaryngol Allied Sci.* 1995;20(6):547-51.
51. Mohide EA, Archibald SD, Tew M, Young JE, Haines T. Postlaryngectomy quality-of-life dimensions identified by patients and health care professionals. *Am J Surg.* 1992;164(6):619-22.
52. Leemans M, van Sluis KE, Van Son R, Van den Brekel MM. Interaction of functional and participation issues on quality of life after total laryngectomy. *Laryngoscope Investig Otolaryngol.* 2020;5(3):453-60.
53. Hilgers FJ, Ackerstaff AH. Comprehensive rehabilitation after total laryngectomy is more than voice alone. *Folia Phoniatr Logop.* 2000;52(1-3):65-73.
54. Brook I. *The Laryngectomee Guide* 2013. 90 p.
55. Zuur JK, Muller SH, de Jongh FH, van Zandwijk N, Hilgers FJ. The physiological rationale of heat and moisture exchangers in post-laryngectomy pulmonary rehabilitation: a review. *Eur Arch Otorhinolaryngol.* 2006;263(1):1-8.
56. Thomachot L, Viviani X, Arnaud S, Boisson C, Martin CD. Comparing two heat and moisture exchangers, one hydrophobic and one hygroscopic, on humidifying efficacy and the rate of nosocomial pneumonia. *Chest.* 1998;114(5):1383-9.
57. Grolman W, Blom ED, Branson RD, Schouwenburg PF, Hamaker RC. An efficiency comparison of four heat and moisture exchangers used in the laryngectomized patient. *Laryngoscope.* 1997;107(6):814-20.
58. Brusasco C, Corradi F, Vargas M, Bona M, Bruno F, Marsili M, et al. In vitro evaluation of heat and moisture exchangers designed for spontaneously breathing tracheostomized patients. *Respir Care.* 2013;58(11):1878-85.
59. Hilgers FJ, Ackerstaff AH, Van As CJ, Balm AJ, Van den Brekel MW, Tan IB. Development and clinical assessment of a heat and moisture exchanger with a multi-magnet automatic tracheostoma valve (Provox FreeHands HME) for vocal and pulmonary rehabilitation after total laryngectomy. *Acta Otolaryngol.* 2003;123(1):91-9.
60. Raynor PC, Kim BG, Ramachandran G, Strommen MR, Horns JH, Streifel AJ. Collection of biological and non-biological particles by new and used filters made from glass and electrostatically charged synthetic fibers. *Indoor Air.* 2008;18(1):51-62.
61. Jones AS, Young PE, Hanafi ZB, Makura ZG, Fenton JE, Hughes JP. A study of the effect of a resistive heat moisture exchanger (Trachinaze) on pulmonary function and blood gas tensions in patients who have undergone a laryngectomy: a randomized control trial of 50 patients studied over a 6-month period. *Head Neck.* 2003;25(5):361-7.
62. Thoman RL, Stoker GL, Ross JC. The efficacy of pursed-lips breathing in patients with chronic obstructive pulmonary disease. *Am Rev Respir Dis.* 1966;93(1):100-6.
63. Tiep BL, Burns M, Kao D, Madison R, Herrera J. Pursed lips breathing training using ear oximetry. *Chest.* 1986;90(2):218-21.

64. Motley HL. The effects of slow deep breathing on the blood gas exchange in emphysema. *Am Rev Respir Dis.* 1963;88:484-92.
65. Mueller RE, Petty TL, Filley GF. Ventilation and arterial blood gas changes induced by pursed lips breathing. *J Appl Physiol.* 1970;28(6):784-9.
66. Schmidt GB, Bombeck CT, Bennett EJ, Kotb KM. Continuous positive airway pressure in the prophylaxis of the adult respiratory distress syndrome (ARDS). *Langenbecks Arch Chir.* 1975;Suppl:439-42.
67. Toremalm NG. A heat-and-moisture exchanger for posttracheotomy care. An experimental study. *Acta Otolaryngol.* 1960;52:461-72.
68. Toremalm NG. The daily amount of tracheo-bronchial secretions in man. A method for continuous tracheal aspiration in laryngectomized and tracheotomized patients. *Acta Otolaryngol Suppl.* 1960;158:43-53.
69. Shanks CA, Sara CA. A reappraisal of the multiple gauze heat and moisture exchanger. *Anaesthesia and intensive care.* 1973;1(5):428-32.
70. Steward DJ. A disposable condenser humidifier for use during anaesthesia. *Can Anaesth Soc J.* 1976;23(2):191-5.
71. Keck T, Durr J, Leiacker R, Rettinger G, Rozsasi A. Tracheal climate in laryngectomees after use of a heat and moisture exchanger. *Laryngoscope.* 2005;115(3):534-7.
72. Vitacca M, Clini E, Foglio K, Scalvini S, Marangoni S, Quadri A, et al. Hygroscopic condenser humidifiers in chronically tracheostomized patients who breathe spontaneously. *Eur Respir J.* 1994;7(11):2026-32.
73. Thomachot L, Viviani X, Arnaud S, Vialet R, Albanese J, Martin C. Preservation of humidity and heat of respiratory gases in spontaneously breathing, tracheostomized patients. *Acta Anaesthesiol Scand.* 1998;42(7):841-4.
74. Nakanishi N, Oto J, Itagaki T, Nakataki E, Onodera M, Nishimura M. Humidification Performance of Passive and Active Humidification Devices Within a Spontaneously Breathing Tracheostomized Cohort. *Respir Care.* 2018;64(2):130-5.
75. Ackerstaff AH, Hilgers FJ, Aaronson NK, Schouwenburg PF, van ZN. [Physical and psychosocial sequelae of total larynx extirpation and the use of a heat and moisture exchanger]. *Ned Tijdschr Geneesk.* 1990;134(50):2438-42.
76. Hilgers FJ, Aaronson NK, Ackerstaff AH, Schouwenburg PF, van Zandwijk N. The influence of a heat and moisture exchanger (HME) on the respiratory symptoms after total laryngectomy. *Clin Otolaryngol Allied Sci.* 1991;16(2):152-6.
77. Ackerstaff AH, Hilgers FJ, Aaronson NK, Balm AJ, Van ZN. Improvements in respiratory and psychosocial functioning following total laryngectomy by the use of a heat and moisture exchanger. *Ann Otol Rhinol Laryngol.* 1993;102(11):878-83.
78. Ackerstaff AH, Hilgers FJ, Aaronson NK, De Boer MF, Meeuwis CA, Knegt PP, et al. Heat and moisture exchangers as a treatment option in the post-operative rehabilitation of laryngectomized patients. *Clin Otolaryngol Allied Sci.* 1995;20(6):504-9.
79. Icuspit P, Yarlagadda B, Garg S, Johnson T, Deschler D. Heat and moisture exchange devices for patients undergoing total laryngectomy. *ORL Head Neck Nurs.* 2014;32(1):20-3.

80. Lorenz KJ, Maier H. [Pulmonary rehabilitation after total laryngectomy using a heat and moisture exchanger (HME)] in German. *Laryngorhinootologie*. 2009;88(8):513-22.
81. Kramp B, Donat M, Dommerich S, Pau HW, Podbielski A. Prospective controlled study of microbial colonization of the trachea in tracheotomized and laryngectomized patients with HME (heat and moisture exchanger). *Acta Otolaryngol*. 2009;129(10):1136-44.
82. Rosso M, Prgomet D, Marjanovic K, Puseljic S, Kraljik N. Pathohistological changes of tracheal epithelium in laryngectomized patients. *Eur Arch Otorhinolaryngol*. 2015;272(11):3539-44.
83. Hilgers FJ, Aaronson NK, Ackerstaff AH, Schouwenburg PF, van ZN. The influence of a heat and moisture exchanger (HME) on the respiratory symptoms after total laryngectomy. *Clin Otolaryngol Allied Sci*. 1991;16(2):152-6.
84. van den Boer C, Muller SH, Vincent AD, Zuchner K, van den Brekel MW, Hilgers FJ. Ex vivo water exchange performance and short-term clinical feasibility assessment of newly developed heat and moisture exchangers for pulmonary rehabilitation after total laryngectomy. *Eur Arch Otorhinolaryngol*. 2014;271(2):359-66.
85. Scheenstra RJ, Muller SH, Vincent A, Ackerstaff AH, Jacobi I, Hilgers FJ. A new heat and moisture exchanger for laryngectomized patients: endotracheal temperature and humidity. *Respir Care*. 2011;56(5):604-11.
86. van den Boer C, Muller SH, Vincent AD, van den Brekel MW, Hilgers FJ. Ex vivo assessment and validation of water exchange performance of 23 heat and moisture exchangers for laryngectomized patients. *Respir Care*. 2014;59(8):1161-71.
87. van den Boer C, Vas Nunes JH, Muller SH, van dN, V, van den Brekel MW, Hilgers FJ. Water Uptake Performance of Hygroscopic Heat and Moisture Exchangers after 24-Hour Tracheostoma Application. *Otolaryngol Head Neck Surg*. 2014;150(6):999-1004.
88. Herranz J, Espino MA, Morado CO. Pulmonary rehabilitation after total laryngectomy: a randomized cross-over clinical trial comparing two different heat and moisture exchangers (HMEs). *Eur Arch Otorhinolaryngol*. 2013;270(9):2479-84.
89. Bien S, Okla S, van As-Brooks CJ, Ackerstaff AH. The effect of a Heat and Moisture Exchanger (Provox HME) on pulmonary protection after total laryngectomy: a randomized controlled study. *Eur Arch Otorhinolaryngol*. 2010;267(3):429-35.
90. Dassonville O, Merol JC, Bozec A, Swierkosz F, Santini J, Chais A, et al. Randomised, multi-centre study of the usefulness of the heat and moisture exchanger (Provox HME(R)) in laryngectomised patients. *Eur Arch Otorhinolaryngol*. 2011;268(11):1647-54.
91. Hennessy M, Bann DV, Patel VA, Saadi R, Krempl GA, Deschler DG, et al. Commentary on the management of total laryngectomy patients during the COVID-19 pandemic. *Head Neck*. 2020.
92. Hess CT. When to use hydrogel dressings. *Adv Skin Wound Care*. 2000;13(1):42.
93. Sarabahi S. Recent advances in topical wound care. *Indian J Plast Surg*. 2012;45(2):379-87.
94. Williams C. Hydrogel dressings: their role in dry and sloughy wounds. *Community Nurse*. 1998;4(9):42-4.

95. Sood A, Granick MS, Tomaselli NL. Wound Dressings and Comparative Effectiveness Data. *Adv Wound Care (New Rochelle)*. 2014;3(8):511-29.
96. Hampton S. A focus on ActiFormCool in the reduction of pain in wounds. *Br J Community Nurs*. 2007;12(9):S37-42.
97. Rippon M, White R, Davies P. Skin adhesives and their role in wound dressings. *Wounds UK*. 2007;3(4):76-86.
98. Blumenstein I, Borger D, Loitsch S, Bott C, Tessmer A, Hartmann F, et al. A glycerin hydrogel-based wound dressing prevents peristomal infections after percutaneous endoscopic gastrostomy (PEG): a prospective, randomized study. *Nutr Clin Pract*. 2012;27(3):422-5.
99. Francesko A, Petkova P, Tzanov T. Hydrogel Dressings for Advanced Wound Management. *Curr Med Chem*. 2017;25(41):5782-97.
100. Scheenstra RJ, Muller SH, Vincent A, Ackerstaff AH, Jacobi I, Hilgers FJ. Short-term endotracheal climate changes and clinical effects of a heat and moisture exchanger with an integrated electrostatic virus and bacterial filter developed for laryngectomized individuals. *Acta Otolaryngol*. 2010;130(6):739-46.
101. Lansaat L, de Kleijn BJ, Hilgers FJ, van der Laan BF, van den Brekel MW. A prospective multicenter clinical feasibility study of a new automatic speaking valve for postlaryngectomy voice rehabilitation. *Eur Arch Otorhinolaryngol*. 2017;274(2):1005-13.
102. Macri GF, Bogaardt H, Parrilla C, Minni A, D'Alatri L, de Vincentiis M, et al. Patients' experiences with HMEs and attachments after total laryngectomy. *Clin Otolaryngol*. 2016;41(6):652-9.
103. Merol JC, Charpiot A, Langagne T, Hemar P, Ackerstaff AH, Hilgers FJ. Randomized controlled trial on postoperative pulmonary humidification after total laryngectomy: External humidifier versus heat and moisture exchanger. *Laryngoscope*. 2012;122(2):275-81.
104. Lorenz KJ, Huverstuhl K, Maier H. [Finger-free speech with the Provox FreeHands HME Automatic Tracheostoma Valve system. Clinical long-term experience]. *HNO*. 2009;57(11):1090-8.
105. Lorenz KJ, Groll K, Ackerstaff AH, Hilgers FJ, Maier H. Hands-free speech after surgical voice rehabilitation with a Provox voice prosthesis: experience with the Provox FreeHands HME tracheostoma valve system. *Eur Arch Otorhinolaryngol*. 2007;264(2):151-7.
106. Ackerstaff AH, Fuller D, Irvin M, Maccracken E, Gaziano J, Stachowiak L. Multicenter study assessing effects of heat and moisture exchanger use on respiratory symptoms and voice quality in laryngectomized individuals. *Otolaryngol Head Neck Surg*. 2003;129(6):705-12.
107. Ward EC, Acton LM, Morton A-L. Stoma Care and Appliances. In: Ward EC, van As-Brooks CJ, editors. *Head and Neck Cancer: Treatment, Rehabilitation, and Outcomes*. San Diego: Plural Publishing; 2007. p. 289-311.
108. Hilgers FJ, Ackerstaff AH. Development and evaluation of a novel tracheostoma button and fixation system (Provox LaryButton and LaryClip adhesive) to facilitate hands-free tracheoesophageal speech. *Acta Otolaryngol*. 2006;126(11):1218-24.
109. Moreno MA, Lewin JS, Hutcheson KA, Bishop Leone JK, Barringer DA, Reece GP. Tracheostomoplasty: A surgical method for improving retention of an intraluminal stoma button for hands-free tracheoesophageal speech. *Head Neck*. 2010;32(12):1674-80.

110. Zuur JK, Muller SH, de Jongh FH, van der Horst MJ, Shehata M, van Leeuwen J, et al. A newly developed tool for intra-tracheal temperature and humidity assessment in laryngectomized individuals: the Airway Climate Explorer (ACE). *Med Biol Eng Comput.* 2007;45(8):737-45.
111. Zuur JK, Muller SH, Vincent A, Sinaasappel M, de Jongh FH, Hilgers FJ. Assessment of tracheal temperature and humidity in laryngectomized individuals and the influence of a heat and moisture exchanger on tracheal climate. *Head Neck.* 2008;30(8):1072-82.
112. Scheenstra RJ, Muller SH, Hilgers FJ. Endotracheal temperature and humidity in laryngectomized patients in a warm and dry environment and the effect of a heat and moisture exchanger. *Head Neck.* 2011;33(9):1285-93.
113. Scheenstra RJ, Muller SH, Vincent A, Sinaasappel M, Zuur JK, Hilgers FJ. Endotracheal temperature and humidity measurements in laryngectomized patients: intra- and inter-patient variability. *Med Biol Eng Comput.* 2009;47(7):773-82.
114. van den Boer C, Muller SH, Vincent AD, Zuchner K, van den Brekel MW, Hilgers FJ. A novel, simplified ex vivo method for measuring water exchange performance of heat and moisture exchangers for tracheostomy application. *Respir Care.* 2013;58(9):1449-58.
115. van den Boer C, Muller SH, van dN, V, Olmos RA, Minni A, Parrilla C, et al. Effects of heat and moisture exchangers on tracheal mucociliary clearance in laryngectomized patients: a multi-center case-control study. *Eur Arch Otorhinolaryngol.* 2015;272(11):3439-50.
116. Hilgers FJ, Ackerstaff AH, Balm AJ, Gregor RT. A new heat and moisture exchanger with speech valve (Provox stomafilter). *Clin Otolaryngol Allied Sci.* 1996;21(5):414-8.
117. Balle VH, Rindso L, Thomsen J. Speech rehabilitation by Provox voice prosthesis combined with heat and moisture exchange filters. *Acta Otolaryngol Suppl.* 1997;529:251-3.
118. Herranz Gonzalez-Botas J, Suarez T, Garcia Carreira B, Martinez Moran A. [Experience with the HME-Provox Stomafilter in laryngectomized patients]. *Acta Otorrinolaringol Esp.* 2001;52(3):221-5.
119. Masson AC, Fouquet ML, Goncalves AJ. Tracheostoma humidifier: influence on secretion and voice of patients with total laryngectomy. *Pro Fono.* 2008;20(3):183-9.
120. Ackerstaff AH, Hilgers FJ, Balm AJ, Tan IB. Long-term compliance of laryngectomized patients with a specialized pulmonary rehabilitation device: Provox Stomafilter. *Laryngoscope.* 1998;108(2):257-60.
121. Scheenstra RJ, Muller SH, Vincent A, Ackerstaff AH, Jacobi I, Hilgers FJ. A new heat and moisture exchanger for laryngectomized patients: endotracheal temperature and humidity. *Respir Care.* 2011;56(5):604-11.
122. Parrilla C, Minni A, Bogaardt H, Macri GF, Battista M, Roukos R, et al. Pulmonary Rehabilitation After Total Laryngectomy: A Multicenter Time-Series Clinical Trial Evaluating the Provox XtraHME in HME-Naive Patients. *Ann Otol Rhinol Laryngol.* 2015;124(9):706-13.
123. Ebersole B, Moran K, Gou J, Ridge J, Schiech L, Liu JC, et al. Heat and moisture exchanger cassettes: Results of a quality/safety initiative to reduce postoperative mucus plugging after total laryngectomy. *Head Neck.* 2020.
124. Ratnayake CBB, Fles R, Tan IB, Baijens LWJ, Pilz W, Meeuwis CA, et al. Multicenter randomized crossover trial evaluating the provox luna in laryngectomized subjects. *Laryngoscope.* 2019;129(10):2354-60.

125. Crawford B, Monz B, Hohlfeld J, Roche N, Rubin B, Magnussen H, et al. Development and validation of a cough and sputum assessment questionnaire. *Respir Med.* 2008;102(11):1545-55.
126. Longobardi Y, Galli J, Di Cesare T, D'Alatri L, Settini S, Mele D, et al. Optimizing Pulmonary Outcomes after Total Laryngectomy: Crossover study on new Heat and Moisture Exchangers. *Otolaryngol Head Neck Surg.* 2022.
127. van As CJ, Hilgers FJ, Koopmans-van Beinum FJ, Ackerstaff AH. The influence of stoma occlusion on aspects of tracheoesophageal voice. *Acta Otolaryngol.* 1998;118(5):732-8.
128. Brook I, Bogaardt H, van As-Brooks C. Long-term use of heat and moisture exchangers among laryngectomees: medical, social, and psychological patterns. *Ann Otol Rhinol Laryngol.* 2013;122(6):358-63.
129. Pedemonte-Sarrias G, Villatoro-Sologaistoa JC, Ale-Inostroza P, Lopez-Vilas M, Leon-Vintro X, Quer-Agusti M. Chronic adherence to heat and moisture exchanger use in laryngectomized patients. *Acta Otorrinolaringol Esp.* 2013;64 4:247-52.
130. Beck ACC, Retel VP, Bunting G, Sethi R, Deschler DG, van den Brekel MWM, et al. Cost-effectiveness analysis of using the heat and moisture exchangers compared with alternative stoma covers in laryngectomy rehabilitation: US perspective. *Head Neck.* 2020.
131. Scheenstra RJ, Muller SH, Vincent A, Ackerstaff AH, Jacobi I, Hilgers FJ. Short-term endotracheal climate changes and clinical effects of a heat and moisture exchanger with an integrated electrostatic virus and bacterial filter developed for laryngectomized individuals. *Acta Otolaryngol.* 2010;130(6):739-46.
132. Kirton OC, DeHaven B, Morgan J, Morejon O, Civetta J. A prospective, randomized comparison of an in-line heat moisture exchange filter and heated wire humidifiers: rates of ventilator-associated early-onset (community-acquired) or late-onset (hospital-acquired) pneumonia and incidence of endotracheal tube occlusion. *Chest.* 1997;112(4):1055-9.
133. Kranabetter R, Leier M, Kammermeier D, Just HM, Heuser D. [The effects of active and passive humidification on ventilation-associated nosocomial pneumonia]. *Anaesthesist.* 2004;53(1):29-35.
134. Subayi L, Chergui K, Beydon L. [Heat and moisture exchanging filters for conditioning of inspired gases in adult anesthesia and resuscitation]. *Ann Fr Anesth Reanim.* 1998;17(7):699-708.
135. Barbara J, Santais MC, Levy DA, Ruff F, Leynadier F. Prevention of latex sensitization in guinea pigs by a bacterial and viral filter used in anaesthesia. *Br J Anaesth.* 2005;95(3):349-54.
136. Boots RJ, Howe S, George N, Harris FM, Faoagali J. Clinical utility of hygroscopic heat and moisture exchangers in intensive care patients. *Crit Care Med.* 1997;25(10):1707-12.
137. Rathgeber J, Kietzmann D, Mergeryan H, Hub R, Zuchner K, Kettler D. Prevention of patient bacterial contamination of anaesthesia-circle-systems: a clinical study of the contamination risk and performance of different heat and moisture exchangers with electret filter (HMEF). *Eur J Anaesthesiol.* 1997;14(4):368-73.
138. Parrinello G, Missale F, Sampieri C, Carobbio ALC, Peretti G. Safe management of laryngectomized patients during the COVID-19 pandemic. *Oral Oncol.* 2020;107:104742.
139. Brook I. Prevention of COVID-19 Infection in Neck Breathers, Including Laryngectomees. *Int Arch Otorhinolaryngol.* 2020;24(3):e253-e4.

140. Yeung DCM, Lai R, Wong EWY, Chan JYK. Care of Patients With a Laryngectomy During the COVID-19 Pandemic. *Otolaryngol Head Neck Surg*. 2020;194599820933185.
141. Lindsley WG, Blachere FM, Beezhold DH, Law BF, Derk RC, Hettick JM, et al. A comparison of performance metrics for cloth face masks as source control devices for simulated cough and exhalation aerosols. *medRxiv*. 2021.
142. Op de Coul BM, Ackerstaff AH, van As-Brooks CJ, van den Hoogen FJ, Meeuwis CA, Manni JJ, et al. Compliance, quality of life and quantitative voice quality aspects of hands-free speech. *Acta Otolaryngol*. 2005;125(6):629-37.
143. Tervonen H, Back L, Juvas A, Rasanen P, Makitie AA, Sintonen H, et al. Automatic speaking valve in speech rehabilitation for laryngectomized patients. *Eur Arch Otorhinolaryngol*. 2005;262(10):816-20.
144. Hamade R, Hewlett N, Scanlon E. A quantitative and qualitative evaluation of an automatic occlusion device for tracheoesophageal speech: the Provox FreeHands HME. *Clin Linguist Phon*. 2006;20(2-3):187-93.
145. Lorenz KJ, Huverstuhl K, Maier H. [Finger-free speech with the Provox FreeHands HME Automatic Tracheostoma Valve system. Clinical long-term experience] in German. *HNO*. 2009;57(11):1090-8.
146. van der Houwen EB, van Kalker TA, Post WJ, Hilgers FJ, van der Laan BF, Verkerke GJ. Does the patch fit the stoma? A study on peristoma geometry and patch use in laryngectomized patients. *Clin Otolaryngol*. 2011;36(3):235-41.
147. Masson AC, Fouquet ML, Goncalves AJ. Tracheostoma humidifier: influence on secretion and voice of patients with total laryngectomy. *Pro Fono*. 2008;20(3):183-9.
148. Lansaat L, Kleijn B, Hilgers F, Laan BV, Brekel MVD. Comparative Study Between Peristomal Patches in Patients with Definitive Tracheostomy. *Int Arch Otorhinolaryngol*. 2018;22(2):130-5.
149. Lewin JS, Lemon J, Bishop-Leone JK, Leyk S, Martin JW, Gillenwater AM. Experience with Barton button and peristomal breathing valve attachments for hands-free tracheoesophageal speech. *Head Neck*. 2000;22(2):142-8.
150. Lewin JS, Montgomery PC, Hutcheson KA, Chambers MS. Further experience with modification of an intraluminal button for hands-free tracheoesophageal speech after laryngectomy. *J Prosthet Dent*. 2009;102(5):328-31.
151. Scheenstra RJ, Muller SH, Hilgers FJ. Endotracheal temperature and humidity in laryngectomized patients in a warm and dry environment and the effect of a heat and moisture exchanger. *Head Neck*. 2011;33(9):1285-93.
152. van den Boer C, Vas Nunes JH, Muller SH, van der Noort V, van den Brekel MW, Hilgers FJ. Water Uptake Performance of Hygroscopic Heat and Moisture Exchangers after 24-Hour Tracheostoma Application. *Otolaryngol Head Neck Surg*. 2014;150(6):999-1004.
153. van den Boer C, Muller SH, van der Noort V, Olmos RA, Minni A, Parrilla C, et al. Effects of heat and moisture exchangers on tracheal mucociliary clearance in laryngectomized patients: a multi-center case-control study. *Eur Arch Otorhinolaryngol*. 2015;272(11):3439-50

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1. Rosso M, Siric L, Kopf T. The role of a heat and moisture exchanger on quality of life after laryngectomy. *Acta Med Croatica*, 75 (2021) 37-40.
2. Searl J, Kearney A, Genoa K, Doyle PC. Clinical Experiences of People With a Laryngectomy During the SARS COVID-19 Pandemic. *Am J Speech Lang Pathol*. 2021:1-16.
3. Ebersole B, Moran K, Gou J, Ridge J, Schiech L, Liu JC, et al. Heat and moisture exchanger cassettes: Results of a quality/safety initiative to reduce postoperative mucus plugging after total laryngectomy. *Head Neck*. 2020.
4. Ratnayake CBB, Fles R, Tan IB, Baijens LWJ, Pilz W, Meeuwis CA, et al. Multicenter randomized crossover trial evaluating the provox Luna® in laryngectomized subjects. *Laryngoscope*. 2019;129(10):2354-60.
5. Lansaat L, Kleijn B, Hilgers F, Laan BV, Brekel MVD. Comparative Study Between Peristomal Patches in Patients with Definitive Tracheostomy. *Int Arch Otorhinolaryngol*. 2018;22(2):130

6. Lansaat L, Boer CV, Muller SH, Noort VV, Brekel MW, Hilgers FJ. Ex vivo humidifying capacity and patient acceptability of stoma cloths in laryngectomized individuals. *Head Neck*. 2017;39(5):921-31.
7. Lansaat L, de Kleijn BJ, Hilgers FJ, van der Laan BF, van den Brekel MW. A prospective multicenter clinical feasibility study of a new automatic speaking valve for postlaryngectomy voice rehabilitation. *Eur Arch Otorhinolaryngol*. 2016;274(2):1005-13.
8. Foreman A, De Santis RJ, Sultanov F, Enepekides DJ, Higgins KM. Heat and moisture exchanger use reduces in-hospital complications following total laryngectomy: a case-control study. *J Otolaryngol Head Neck Surg*. 2016;45(1):40.
9. Macri GF, Bogaardt H, Parrilla C, Minni A, D'Alatri L, de Vincentiis M, et al. Patients' experiences with HMEs and attachments after total laryngectomy. *Clin Otolaryngol*. 2016;41(6):652-9.
10. Parrilla C, Minni A, Bogaardt H, Macri GF, Battista M, Roukos R, et al. Pulmonary Rehabilitation After Total Laryngectomy: A Multicenter Time-Series Clinical Trial Evaluating the Provox XtraHME in HME-Naive Patients. *Ann Otol Rhinol Laryngol*. 2015;124(9):706-13.
11. Rosso M, Prgomet D, Marjanovic K, Puseljic S, Kraljik N. Pathohistological changes of tracheal epithelium in laryngectomized patients. *Eur Arch Otorhinolaryngol*. 2015;272(11):3539-44.
12. van den Boer C, Muller SH, van dN, V, Olmos RA, Minni A, Parrilla C, et al. Effects of heat and moisture exchangers on tracheal mucociliary clearance in laryngectomized patients: a multi-center case-control study. *Eur Arch Otorhinolaryngol*. 2015;272(11):3439-50.
13. van den Boer C, van Harten MC, Hilgers FJ, van den Brekel MW, Retel VP. Incidence of severe tracheobronchitis and pneumonia in laryngectomized patients: a retrospective clinical study and a European-wide survey among head and neck surgeons. *Eur Arch Otorhinolaryngol*. 2014;271(12):3297-303.
14. Brook I, Bogaardt H, van As-Brooks C. Long-term use of heat and moisture exchangers among laryngectomees: medical, social, and psychological patterns. *Ann Otol Rhinol Laryngol*. 2013;122(6):358-63.
15. Herranz J, Espino MA, Morado CO. Pulmonary rehabilitation after total laryngectomy: a randomized cross-over clinical trial comparing two different heat and moisture exchangers (HMEs). *Eur Arch Otorhinolaryngol*. 2013;270(9):2479-84.
16. Dirven R, Clark JR, Wismans JG, McGuinness J, Palme CE, Blyth K, et al. A new custom moldable external neck brace (ENB 2.0) to improve hands-free speech in laryngectomized patients. *Laryngoscope*. 2013;123(9):2209-15.
17. Pedemonte-Sarrias G, Villatoro-Sologaistoa JC, Ale-Inostroza P, Lopez-Vilas M, Leon-Vintro X, Quer-Agusti M. Chronic adherence to heat and moisture exchanger use in laryngectomized patients. *Acta Otorrinolaringol Esp*. 2013;64 4:247-52.
18. Hilgers FJ, Dirven R, Wouters Y, Jacobi I, Marres HA, van den Brekel MW. A multicenter, prospective, clinical trial evaluating a novel adhesive baseplate (Provox StabiliBase™) for peristomal attachment of postlaryngectomy pulmonary and voice rehabilitation devices. *Laryngoscope*. 2012;122(11):2447-53.

19. Dirven R, Kooijman PG, Wouters Y, Marres HA. Clinical use of a neck brace to improve hands-free speech in laryngectomized patients. *Laryngoscope*. 2012.
20. Merol JC, Charpiot A, Langagne T, Hemar P, Ackerstaff AH, Hilgers FJ. Randomized controlled trial on postoperative pulmonary humidification after total laryngectomy: External humidifier versus heat and moisture exchanger. *Laryngoscope*. 2012;122(2):275-81.
21. Scheenstra RJ, Muller SH, Hilgers FJ. Endotracheal temperature and humidity in laryngectomized patients in a warm and dry environment and the effect of a heat and moisture exchanger. *Head Neck*. 2011;33(9):1285-93.
22. Dassonville O, Merol JC, Bozec A, Swierkosz F, Santini J, Chais A, et al. Randomised, multi-centre study of the usefulness of the heat and moisture exchanger (Provox HME(R)) in laryngectomised patients. *Eur Arch Otorhinolaryngol*. 2011;268(11):1647-54.
23. Scheenstra RJ, Muller SH, Vincent A, Ackerstaff AH, Jacobi I, Hilgers FJ. A new heat and moisture exchanger for laryngectomized patients: endotracheal temperature and humidity. *Respir Care*. 2011;56(5):604-11.
24. Scheenstra RJ, Muller SH, Vincent A, Ackerstaff AH, Jacobi I, Hilgers FJ. Short-term endotracheal climate changes and clinical effects of a heat and moisture exchanger with an integrated electrostatic virus and bacterial filter developed for laryngectomized individuals. *Acta Otolaryngol*. 2010;130(6):739-46.
25. Scheenstra RJ, Muller SH, Vincent A, Sinaasappel M, Hilgers FJ. Influence of breathing resistance of heat and moisture exchangers on tracheal climate and breathing pattern in laryngectomized individuals. *Head Neck*. 2010;32(8):1069-78.
26. Lorenz KJ, Huverstuhl K, Maier H. [Finger-free speech with the Provox FreeHands HME Automatic Tracheostoma Valve system. Clinical long-term experience] in German. *HNO*. 2009;57(11):1090-8.
27. Bien S, Okla S, van As-Brooks CJ, Ackerstaff AH. The effect of a Heat and Moisture Exchanger (Provox HME) on pulmonary protection after total laryngectomy: a randomized controlled study. *Eur Arch Otorhinolaryngol*. 2010;267(3):429-35.
28. Zuur JK, Muller SH, Vincent A, Sinaasappel M, de Jongh FH, Hilgers FJ. The influence of a heat and moisture exchanger on tracheal climate in a cold environment. *Med Eng Phys*. 2009;31(7):852-7.
29. Scheenstra RJ, Muller SH, Vincent A, Sinaasappel M, Zuur JK, Hilgers FJ. Endotracheal temperature and humidity measurements in laryngectomized patients: intra- and inter-patient variability. *Med Biol Eng Comput*. 2009;47(7):773-82.
30. Kramp B, Donat M, Dommerich S, Pau HW, Podbielski A. Prospective controlled study of microbial colonization of the trachea in tracheotomized and laryngectomized patients with HME (heat and moisture exchanger). *Acta Otolaryngol*. 2009;129(10):1136-44.
31. Masson AC, Fouquet ML, Goncalves AJ. Tracheostoma humidifier: influence on secretion and voice of patients with total laryngectomy. *Pro Fono*. 2008;20(3):183-9.

32. Zuur JK, Muller SH, Vincent A, Sinaasappel M, de Jongh FH, Hilgers FJ. Assessment of tracheal temperature and humidity in laryngectomized individuals and the influence of a heat and moisture exchanger on tracheal climate. *Head Neck*. 2008;30(8):1072-82.
33. Dupuis P, Guertin L, Rainville MS, Prud'homme DL, Lavigne F. Montreal's experience with Cyranose heat and moisture exchanger use in 15 laryngectomized patients. *J Otolaryngol*. 2007;36(4):208-12.
34. Zuur JK, Muller SH, Sinaasappel M, Hart GA, Van ZN, Hilgers FJ. Influence of heat and moisture exchanger respiratory load on transcutaneous oxygenation in laryngectomized individuals: a randomized crossover study. *Head Neck*. 2007;29(12):1102-10.
35. Hilgers FJ, Ackerstaff AH. Development and evaluation of a novel tracheostoma button and fixation system (Provox LaryButton and LaryClip adhesive) to facilitate hands-free tracheoesophageal speech. *Acta Otolaryngol*. 2006;126(11):1218-24.
36. Lorenz KJ, Groll K, Ackerstaff AH, Hilgers FJ, Maier H. Hands-free speech after surgical voice rehabilitation with a Provox voice prosthesis: experience with the Provox FreeHands HME tracheostoma valve system. *Eur Arch Otorhinolaryngol*. 2007;264(2):151-7.
37. Hamade R, Hewlett N, Scanlon E. A quantitative and qualitative evaluation of an automatic occlusion device for tracheoesophageal speech: the Provox FreeHands HME. *Clin Linguist Phon*. 2006;20(2-3):187-93.
38. Op de Coul BM, Ackerstaff AH, van As-Brooks CJ, van den Hoogen FJ, Meeuwis CA, Manni JJ, et al. Compliance, quality of life and quantitative voice quality aspects of hands-free speech. *Acta Otolaryngol*. 2005;125(6):629-37.
39. Keck T, Durr J, Leiacker R, Rettinger G, Rozsasi A. Tracheal climate in laryngectomees after use of a heat and moisture exchanger. *Laryngoscope*. 2005;115(3):534-7.
40. Tervonen H, Back L, Juvas A, Rasanen P, Makitie AA, Sintonen H, et al. Automatic speaking valve in speech rehabilitation for laryngectomized patients. *Eur Arch Otorhinolaryngol*. 2005;262(10):816-20.
41. Ackerstaff AH, Fuller D, Irvin M, Maccracken E, Gaziano J, Stachowiak L. Multicenter study assessing effects of heat and moisture exchanger use on respiratory symptoms and voice quality in laryngectomized individuals. *Otolaryngol Head Neck Surg*. 2003;129(6):705-12.
42. Moerman M, Lawson G, Andry G, Remacle M. The Belgian experience with the cyranose heat moisture exchange filter. A multicentric pilot study of 12 total laryngectomees. *Eur Arch Otorhinolaryngol*. 2003;260(6):301-3.
43. Hilgers FJ, Ackerstaff AH, Van As CJ, Balm AJ, Van den Brekel MW, Tan IB. Development and clinical assessment of a heat and moisture exchanger with a multi-magnet automatic tracheostoma valve (Provox FreeHands HME) for vocal and pulmonary rehabilitation after total laryngectomy. *Acta Otolaryngol*. 2003;123(1):91-9.
44. Herranz Gonzalez-Botas J, Suarez T, Garcia CB, Martinez MA. [Experience with the HME-Provox Stomafilter in laryngectomized patients] in Spanish. *Acta Otorrinolaringol Esp*. 2001;52(3):221-5.

45. Van As CJ, Hilgers FJ, Koopmans-van Beinum FJ, Ackerstaff AH. The influence of stoma occlusion on aspects of tracheoesophageal voice. *Acta Otolaryngol.* 1998;118(5):732-8.
46. Ackerstaff AH, Hilgers FJ, Balm AJ, Tan IB. Long-term compliance of laryngectomized patients with a specialized pulmonary rehabilitation device: Provox Stomafilter. *Laryngoscope.* 1998;108(2):257-60.
47. Hilgers FJM, Ackerstaff AH, Balm AJ, Gregor RT. A new heat and moisture exchanger with speech valve (Provox stomafilter). *Clin Otolaryngol Allied Sci.* 1996;21(5):414-8.
48. Ackerstaff AH, Hilgers FJ, Aaronson NK, De Boer MF, Meeuwis CA, Knegt PP, et al. Heat and moisture exchangers as a treatment option in the post-operative rehabilitation of laryngectomized patients. *Clin Otolaryngol Allied Sci.* 1995;20(6):504-9.
49. Ackerstaff AH, Hilgers FJ, Aaronson NK, Balm AJ, Van ZN. Improvements in respiratory and psychosocial functioning following total laryngectomy by the use of a heat and moisture exchanger. *Ann Otol Rhinol Laryngol.* 1993;102(11):878-83.
50. Hilgers FJ, Aaronson NK, Ackerstaff AH, Schouwenburg PF, van ZN. The influence of a heat and moisture exchanger (HME) on the respiratory symptoms after total laryngectomy. *Clin Otolaryngol Allied Sci.* 1991;16(2):152-6.
51. Ackerstaff AH, Hilgers FJ, Aaronson NK, Schouwenburg PF, Van ZN. [Physical and psychosocial sequelae of total larynx extirpation and the use of a heat and moisture exchanger] in Dutch. *Ned Tijdschr Geneesk.* 1990;134(50):2438-42.
52. McRae D, Young P, Hamilton J, Jones A. Raising airway resistance in laryngectomees increases tissue oxygen saturation. *Clin Otolaryngol Allied Sci.* 1996;21(4):366-8.
53. Jones AS, Young PE, Hanafi ZB, Makura ZG, Fenton JE, Hughes JP. A study of the effect of a resistive heat moisture exchanger (Trachinaze) on pulmonary function and blood gas tensions in patients who have undergone a laryngectomy: a randomized control trial of 50 patients studied over a 6-month period. *Head Neck.* 2003;25(5):361-7.

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Atos Medical AB, Hyllie Boulevard 17, SE-215 32 Malmö, Suécia

Fabricante:

Atos Medical AB, Kraftgatan 8, SE-242 35 Hörby, Suécia Tel: +46 (0) 415 198 00

Email: info@atosmedical.com