
# TheraBite

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# TheraBite®

## Jaw Motion Rehabilitation System™

## Sumário

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1. van der Molen L, van Rossum MA, Burkhead LM, Smeele LE, Rasch CR, Hilgers FJ. A randomized preventive rehabilitation trial in advanced head and neck cancer patients treated with chemoradiotherapy: feasibility, compliance, and short-term effects. *Dysphagia*. 2011;26(2):155-70.
2. Dijkstra PU, Kalk WW, Roodenburg JL. Trismus in head and neck oncology: a systematic review. *Oral Oncol*. 2004;40(9):879-89.
3. Dijkstra PU, Huisman PM, Roodenburg JL. Criteria for trismus in head and neck oncology. *Int J Oral Maxillofac Surg*. 2006;35(4):337-42.
4. Dhanrajani PJ, Jonaidel O. Trismus: aetiology, differential diagnosis and treatment. *Dent Update*. 2002;29(2):88-92, 4.
5. Scott B, Butterworth C, Lowe D, Rogers SN. Factors associated with restricted mouth opening and its relationship to health-related quality of life in patients attending a Maxillofacial Oncology clinic. *Oral Oncol*. 2008;44(5):430-8.
6. Ozyar E, Cengiz M, Gurkaynak M, Atahan IL. Trismus as a presenting symptom in nasopharyngeal carcinoma. *Radiother Oncol*. 2005;77(1):73-6.
7. Cohen EG, Deschler DG, Walsh K, Hayden RE. Early use of a mechanical stretching device to improve mandibular mobility after composite resection: a pilot study. *Arch Phys Med Rehabil*. 2005;86(7):1416-9.
8. Ichimura K, Tanaka T. Trismus in patients with malignant tumours in the head and neck. *J Laryngol Otol*. 1993;107(11):1017-20.
9. Nguyen TD, Panis X, Froissart D, Legros M, Coninx P, Loirette M. Analysis of late complications after rapid hyperfractionated radiotherapy in advanced head and neck cancers. *Int J Radiat Oncol Biol Phys*. 1988;14(1):23-5.
10. van der Molen L, vanRossum M, Smeele L, Rasch C, Hilgers F. A randomized controlled trial on preventive swallowing exercises in advanced head and neck cancer treated with chemoradiotherapy: 1-year functional outcomes. Preventive rehabilitation in patients treated with chemoradiation for advanced head and neck cancer. Thesis ed. Amsterdam: Nederlands Kanker Instituut; 2011.

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11. van der Geer SJ, van Rijn PV, Kamstra JI, Roodenburg JLN, Dijkstra PU. Criterion for trismus in head and neck cancer patients: a verification study. *Support Care Cancer*. 2018;27(3):1129-37.
12. Beekhuis G, Harrington E. Trismus. Etiology and management of inability to open the mouth. *Laryngoscope*. 1965;75:1234-58.
13. Kent LM, Brennan MT, Noll JL, Fox PC, Burri SH, Hunter JC, et al. Radiation-induced trismus in head and neck cancer patients. *Support Care Cancer*. 2008;16(3):305-9.
14. Johnson J, van As-Brooks CJ, Fagerberg-Mohlin B, Finizia C. Trismus in head and neck cancer patients in Sweden: incidence and risk factors. *Med Sci Monit*. 2010;16(6):CR278-CR82.
15. Goldstein M, Maxymiw WG, Cummings BJ, Wood RE. The effects of antitumor irradiation on mandibular opening and mobility: a prospective study of 58 patients. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod*. 1999;88(3):365-73.
16. Horst RW. Trismus: its causes, effects and treatment. *ORL Head Neck Nurs*. 1994;12(2):11-2.
17. Epstein JB, Emerton S, Kolbinson DA, Le ND, Phillips N, Stevenson-Moore P, et al. Quality of life and oral function following radiotherapy for head and neck cancer. *Head Neck*. 1999;21(1):1-11.
18. Weber C, Dommerich S, Pau HW, Kramp B. Limited mouth opening after primary therapy of head and neck cancer. *Oral Maxillofac Surg*. 2010;14(3):169-73.
19. Nordgren M, Abendstein H, Jannert M, Boysen M, Ahlner-Elmqvist M, Silander E, et al. Health-related quality of life five years after diagnosis of laryngeal carcinoma. *Int J Radiat Oncol Biol Phys*. 2003;56(5):1333-43.
20. Hammerlid E, Silander E, Hornestam L, Sullivan M. Health-related quality of life three years after diagnosis of head and neck cancer--a longitudinal study. *Head Neck*. 2001;23(2):113-25.
21. Lee R, Slevin N, Musgrove B, Swindell R, Molassiotis A. Prediction of post-treatment trismus in head and neck cancer patients. *Br J Oral Maxillofac Surg*. 2012;50(4):328-32.
22. Lee LY, Chen SC, Chen WC, Huang BS, Lin CY. Postradiation trismus and its impact on quality of life in patients with head and neck cancer. *Oral Surg Oral Med Oral Pathol Oral Radiol*. 2015;119(2):187-95.

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23. Kondo T, Sugauchi A, Yabuno Y, Kobashi H, Amano K, Aikawa T, et al. Performance status scale for head and neck scores for oral cancer survivors: predictors and factors for improving quality of life. *Clin Oral Investig*. 2019;23(4):1575-82.
24. Malkawi Z, Al-Omiri MK, Khraisat A. Risk indicators of postoperative complications following surgical extraction of lower third molars. *Med Princ Pract*. 2011;20(4):321-5.
25. van der Geer SJ, van Rijn PV, Kamstra JI, Langendijk JA, van der Laan B, Roodenburg JLN, et al. Prevalence and prediction of trismus in patients with head and neck cancer: A cross-sectional study. *Head Neck*. 2019;41(1):64-71.
26. Rapidis AD, Dijkstra PU, Roodenburg JL, Rodrigo JP, Rinaldo A, Strojan P, et al. Trismus in patients with head and neck cancer: etiopathogenesis, diagnosis and management. *Clin Otolaryngol*. 2015;40(6):516-26.
27. McCaul LK. Oral and dental management for head and neck cancer patients treated by chemotherapy and radiotherapy. *Dent Update*. 2012;39(2):135-8, 40.
28. Weber C. Investigations into maximal interincisor opening in patients with malignant oropharyngeal tumours after completion of primary therapy. Oral Presentation; German National ENT Congress; 2006.
29. Whitmyer CC, Waskowski JC, Iffland HA. Radiotherapy and oral sequelae: preventive and management protocols. *J Dent Hyg*. 1997;71(1):23-9.
30. Hashibe M, Ritz B, Le AD, Li G, Sankaranarayanan R, Zhang ZF. Radiotherapy for oral cancer as a risk factor for second primary cancers. *Cancer Lett*. 2005;220(2):185-95.
31. Jeremic G, Venkatesan V, Hallock A, Scott D, Hammond A, Read N, et al. Trismus following treatment of head and neck cancer. *J Otolaryngol Head Neck Surg*. 2011;40(4):323-9.
32. van der Geer SJ, Kamstra JI, Roodenburg JL, van Leeuwen M, Reintsema H, Langendijk JA, et al. Predictors for trismus in patients receiving radiotherapy. *Acta Oncol*. 2016;55(11):1318-23.
33. Sebastian MH, Moffett BC. The effects of continuous passive motion on the temporomandibular joint after surgery. Part II. Appliance improvement, normal subject evaluation, pilot clinical trial. *Oral Surg Oral Med Oral Pathol*. 1989;67(6):644-53.
34. Standley PR, Meltzer K. In vitro modeling of repetitive motion strain and manual medicine treatments: potential roles for pro- and anti-inflammatory cytokines. *J Bodyw Mov Ther*. 2008;12(3):201-3.

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35. Kamstra JI, Roodenburg JL, Beurskens CH, Reintsema H, Dijkstra PU. TheraBite exercises to treat trismus secondary to head and neck cancer. *Support Care Cancer*. 2013;21(4):951-7.
36. Bensadoun RJ, Riesenbeck D, Lockhart PB, Elting LS, Spijkervet FK, Brennan MT. A systematic review of trismus induced by cancer therapies in head and neck cancer patients. *Support Care Cancer*. 2010;18(8):1033-8.
37. Ezzat YE, Sharka RM, Huzaimi AA, Al-Zahrani KM, Abed HH. The role of exercise therapy in managing post-radiotherapy trismus in head and neck cancer. *J Taibah Univ Med Sci*. 2021;16(1):127-33.
38. Melchers LJ, Van WE, Beurskens CH, Reintsema H, Slagter AP, Roodenburg JL, et al. Exercise adherence in patients with trismus due to head and neck oncology: a qualitative study into the use of the Therabite((R)). *Int J Oral Maxillofac Surg*. 2009;38(9):947-54.
39. Dijkstra PU, Sterken MW, Pater R, Spijkervet FK, Roodenburg JL. Exercise therapy for trismus in head and neck cancer. *Oral Oncol*. 2007;43(4):389-94.
40. Buchbinder D, Currivan RB, Kaplan AJ, Urken ML. Mobilization regimens for the prevention of jaw hypomobility in the radiated patient: a comparison of three techniques. *J Oral Maxillofac Surg*. 1993;51(8):863-7.
41. Pauli N, Fagerberg-Mohlin B, Andrell P, Finizia C. Exercise intervention for the treatment of trismus in head and neck cancer. *Acta Oncol*. 2014;53(4):502-9.
42. Lee R, Yeo ST, Rogers SN, Caress AL, Molassiotis A, Ryder D, et al. Randomised feasibility study to compare the use of Therabite((R)) with wooden spatulas to relieve and prevent trismus in patients with cancer of the head and neck. *Br J Oral Maxillofac Surg*. 2018;56(4):283-91.
43. van der Geer SJ, Reintsema H, Kamstra JI, Roodenburg JLN, Dijkstra PU. The use of stretching devices for treatment of trismus in head and neck cancer patients: a randomized controlled trial. *Support Care Cancer*. 2020;28(1):9-11.
44. Messing BP, Saunders JR, Dietrich-Burns K, Farrell SS, Pyrke-Fairchild J, Melick CF, et al. Trismus management in the head and neck cancer patient undergoing organ preservation treatment. 2009.
45. Montalvo C, Finizia C, Pauli N, Fagerberg-Mohlin B, Andrell P. Impact of exercise with TheraBite device on trismus and health-related quality of life: A prospective study. *Ear Nose Throat J*. 2017;96(1):E1-E6.

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46. Pauli N, Svensson U, Karlsson T, Finizia C. Exercise intervention for the treatment of trismus in head and neck cancer - a prospective two-year follow-up study. *Acta Oncol.* 2016;55(6):1-7.
47. Chee S, Byrnes YM, Chorath KT, Rajasekaran K, Deng J. Interventions for Trismus in Head and Neck Cancer Patients: A Systematic Review of Randomized Controlled Trials. *Integr Cancer Ther.* 2021;20:15347354211006474.
48. Scherpenhuizen A, van Waes AM, Janssen LM, Van Cann EM, Stegeman I. The effect of exercise therapy in head and neck cancer patients in the treatment of radiotherapy- induced trismus: A systematic review. *Oral Oncol.* 2015;51(8):745-50.
49. Karlis V, Glickman RS. Evaluation of mobilization regimens for patients with trismus. Poster presentation; Annual meeting of the American Association of Oral And Maxillofacial Surgeons; 2001.
50. Melchers LJ, Van WE, Beurskens CH, Reintsema H, Slagter AP, Roodenburg JL, et al. Exercise adherence in patients with trismus due to head and neck oncology: a qualitative study into the use of the Therabite. *Int J Oral Maxillofac Surg.* 2009;38(9):947-54.
51. Bensadoun RJ, Riesenbeck D, Lockhart PB, Elting LS, Spijkervet FK, Brennan MT. A systematic review of trismus induced by cancer therapies in head and neck cancer patients. *Support Care Cancer.* 2010;18(8):1033-8.
52. Tang Y, Shen Q, Wang Y, Lu K, Wang Y, Peng Y. A randomized prospective study of rehabilitation therapy in the treatment of radiation-induced dysphagia and trismus. *Strahlenther Onkol.* 2011;187(1):39-44.
53. Treister N, Duncan C, Cutler C, Lehmann L. How we treat oral chronic graft-versus- host disease. *Blood.* 2012;120(17):3407-18.
54. Pauli N, Andrell P, Johansson M, Fagerberg-Mohlin B, Finizia C. Treating trismus: A prospective study on effect and compliance to jaw exercise therapy in head and neck cancer. *Head Neck.* 2015;37(12):1738-44.
55. Lee R, Yeo ST, Rogers SN, Caress AL, Molassiotis A, Ryder D, et al. Randomised feasibility study to compare the use of Therabite((R)) with wooden spatulas to relieve and prevent trismus in patients with cancer of the head and neck. *The British journal of oral & maxillofacial surgery.* 2018 May;56(4):283-291.
56. Montalvo C, Finizia C, Pauli N, Fagerberg-Mohlin B, Andrell P. Impact of exercise with TheraBite device on trismus and health-related quality of life: A prospective study. *Ear Nose Throat J.* 2020:145561320961727.

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57. Zatarain LA, Smith DK, Deng J, Gilbert J, Dietrich MS, Niemann KJ, et al. A Randomized Feasibility Trial to Evaluate Use of the Jaw Dynasplint to Prevent Trismus in Patients With Head and Neck Cancer Receiving Primary or Adjuvant Radiation-Based Therapy. *Integr Cancer Ther.* 2018;1534735418784363.

58. Bragante KC, Groisman S, Carboni C, Baiocchi JMT, da Motta NW, Silva MF, et al. Efficacy of exercise therapy during radiotherapy to prevent reduction in mouth opening in patients with head and neck cancer: A randomized controlled trial. *Oral Surg Oral Med Oral Pathol Oral Radiol.* 2020;129(1):27-38.

59. Loorents V, Rosell J, Karlsson C, Lidback M, Hultman K, Borjeson S. Prophylactic training for the prevention of radiotherapy-induced trismus - a randomised study. *Acta Oncol.* 2014;53(4):530-8.

60. Gassner RJ, Buckley MJ, Studer RK, Evans CH, Agarwal S. Interaction of strain and interleukin-1 in articular cartilage: effects on proteoglycan synthesis in chondrocytes. *Int J Oral Maxillofac Surg.* 2000;29(5):389-94.

61. Robbins JL. Rehabilitation after temporomandibular joint surgery: a review of the literature and guidelines for practice. *TMJournal.* 2002;2(1):39-41.

62. Maloney GE, Mehta N, Forgione AG, Zawawi KH, Al-Badawi EA, Driscoll SE. Effect of a passive jaw motion device on pain and range of motion in TMD patients not responding to flat plane intraoral appliances. *Cranio.* 2002;20(1):55-66.

63. Karlis V, Glickman RS. Evaluation of efficacy of regimens for closed lock patients. Oral presentation; Annual meeting of the American Association of Oral And Maxillofacial Surgeons; 1994.

64. Heres Diddens A, Kraaijenga S, Coupe V, Hilgers F, van der Molen L, Smeele L, et al. The cost-effectiveness of TheraBite(R) as treatment for acute myogenic temporomandibular disorder. *Cranio.* 2016;35(5):1-8.

65. Reed DN, Hall DL, Cottle JH, Frimenko K, Horton CK, Abu Sharkh F, et al. Dental management of scleroderma patients using pentoxifylline plus vitamin E with and without TheraBite((R)) to reduce trismus: Two case reports and brief review of literature. *Clin Case Rep.* 2020;8(2):247-53.

66. Yim CR, Uhrich RM, Perez L, Jr. Treatment of Temporomandibular Joint Heterotopic Ossification: A Novel Protocol With Multimodal Therapy Based on Literature Review and Presentation of a Unique Case Report. *Journal of oral and maxillofacial surgery : official journal of the American Association of Oral and Maxillofacial Surgeons.* 2022;80(5):869-88.

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67. Salter RB. The biologic concept of continuous passive motion of synovial joints. The first 18 years of basic research and its clinical application. *Clin Orthop Relat Res.* 1989(242):12-25.
68. Salter RB. The physiologic basis of continuous passive motion for articular cartilage healing and regeneration. *Hand Clin.* 1994;10(2):211-9.
69. Kim HK, Kerr RG, Cruz TF, Salter RB. Effects of continuous passive motion and immobilization on synovitis and cartilage degradation in antigen induced arthritis. *J Rheumatol.* 1995;22(9):1714-21.
70. Salter RB. History of rest and motion and the scientific basis for early continuous passive motion. *Hand Clin.* 1996;12(1):1-11.
71. Nicolakis P, Erdogmus B, Kopf A, Djaber-Ansari A, Piehslinger E, Fialka-Moser V. Exercise therapy for craniomandibular disorders. *Arch Phys Med Rehabil.* 2000;81(9):1137- 42.
72. Gassner RJ, Agarwal S. Drs. Gassner and Agarwal respond to Maloney, et al.'s article in the January 2002 issue of CRANIO. *Cranio.* 2002;20(3):152-3.
73. Kraaijenga S, van der Molen L, van TH, Hilgers F, Smeele L. Treatment of myogenic temporomandibular disorder: a prospective randomized clinical trial, comparing a mechanical stretching device (TheraBite(R)) with standard physical therapy exercise. *Cranio.* 2014;32(3):208-16.
74. Lloyd Morris EH, Estilow T, Glanzman AM, Cusack SV, Yum SW. Improving Temporomandibular Range of Motion in People With Duchenne Muscular Dystrophy and Spinal Muscular Atrophy. *Am J Occup Ther.* 2020;74(2):7402205080p1-p10.
75. Sydow E, Van der Elst K, Verschueren P, Lenaerts J, Westhovens R, De Langhe E. A mixed method study exploring patient satisfaction and feasibility of two exercise programmes in systemic sclerosis-associated microstomia. *Rheumatol Adv Pract.* 2021;5(2):rkab017.
76. Kaban LB, Perrott DH, Fisher K. A protocol for management of temporomandibular joint ankylosis. *J Oral Maxillofac Surg.* 1990;48(11):1145-51.
77. Stack BC, Gregory E, Gjerde G, Hanssen JI, Leivseth G. Refractive patients: TMJ splint therapy combined with open joint surgery: 60 patients; 117 joints. *Funct Orthod.* 2001;18(4):24-36.
78. Warwas FB, Wiedemeyer V, Kramer FJ, Heim N. Trismus Due to Hyperplastic Coronoid Process: Series of Five Cases and Review of Literature. *The Journal of craniofacial surgery.* 2021;32(8):2798-801.

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79. Aryanpour Z, Maglic D, Zahid R, Tuncer FB, Gociman BR, Siddiqi FA. Mandibular Myofibroma and Severe Trismus: A Complex Case and Review of Complications. *Plast Reconstr Surg Glob Open*. 2022;10(6):e4380.

80. Carnaby-Mann G, Crary MA, Schmalfluss I, Amdur R. "Pharyngocise": randomized controlled trial of preventative exercises to maintain muscle structure and swallowing function during head-and-neck chemoradiotherapy. *Int J Radiat Oncol Biol Phys*. 2012;83(1):210-9.

81. Burkhead LM, Sapienza CM, Rosenbek JC. Strength-training exercise in dysphagia rehabilitation: principles, procedures, and directions for future research. *Dysphagia*. 2007;22(3):251-65.

82. Burkhead L, Rosenbek J, Sapienza C, Musson N, Davenport P, Schechtman O. The effect of jaw and tongue position on suprahyoid muscle function during swallowing. *Dysphagia*. 2006;21(4):298.

83. Burkhead LM, Rosenbek JC, Sapienza CM, Musson N, Davenport P, Schechtman O. The effect of Jaw and Tongue position on suprahyoid muscle functioning during swallowing. Oral presentation; EGDG congress, Malmo; 2006.

84. Karsten RT, Chargi N, van der Molen L, van Son R, de Bree R, Al-Mamgani A, et al. Dysphagia, trismus and speech impairment following radiation-based treatment for advanced stage oropharyngeal carcinoma: a one-year prospective evaluation. *Eur Arch Otorhinolaryngol*. 2022;279(2):1003-27.

85. Wang YH, Huang YA, Chen IH, Hou WH, Kang YN. Exercise for Trismus Prevention in Patients with Head and Neck Cancer: A Network Meta-Analysis of Randomized Controlled Trials. *Healthcare (Basel)*. 2022;10(3).

86. Shao CH, Chiang CC, Huang TW. Exercise therapy for cancer treatment-induced trismus in patients with head and neck cancer: A systematic review and meta-analysis of randomized controlled trials. *Radiother Oncol*. 2020;151:249-55.

87. Kamstra JI, van Leeuwen M, Roodenburg JL, Dijkstra PU. Exercise therapy for trismus secondary to head and neck cancer: A systematic review. *Head Neck*. 2016;39(1):160-9.

88. van der Molen L, van Rossum MA, Rasch CR, Smeele LE, Hilgers FJ. Two-year results of a prospective preventive swallowing rehabilitation trial in patients treated with chemoradiation for advanced head and neck cancer. *Eur Arch Otorhinolaryngol*. 2014;271(5):1257-70.

89. Karsten RT, van der Molen L, Hamming-Vrieze O, van Son R, Hilgers FJM, van den Brekel MWM, et al. Long-term swallowing, trismus, and speech outcomes after combined

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chemoradiotherapy and preventive rehabilitation for head and neck cancer; 10-year plus update. *Head Neck*. 2020 Aug;42(8):1907-1918.

90. Kraaijenga SA, van der Molen L, Jacobi I, Hamming-Vrieze O, Hilgers FJ, van den Brekel MW. Prospective clinical study on long-term swallowing function and voice quality in advanced head and neck cancer patients treated with concurrent chemoradiotherapy and preventive swallowing exercises. *Eur Arch Otorhinolaryngol*. 2015;272(11):3521-31.

91. Messing BP, Ryniack K, Shellenberg J, Hieber B, Ulmer K, Gold D, et al. Multidisciplinary Care of the Head and Neck Cancer Patient. *Perspectives on Voice and Voice Disorders*. 2012;22:14-24.

92. Retel VP, van der Molen L, Steuten LM, van den Brekel MW, Hilgers FJ. A costeffectiveness analysis of using TheraBite in a preventive exercise program for patients with advanced head and neck cancer treated with concomitant chemo-radiotherapy. *Eur Arch Otorhinolaryngol*. 2016;273(3):709-18.

93. Karlsson O, Karlsson T, Pauli N, Andrell P, Finizia C. Jaw exercise therapy for the treatment of trismus in head and neck Cancer: a prospective three-year follow-up study. *Support Care Cancer*. 2021 Jul;29(7):3793-3800.

94. E C, M D, K C, V A, P M, C F, et al. Trismus therapy devices: A systematic review. *Oral Oncol*. 2022;126:105728.

95. Eisbruch A, Schwartz M, Rasch C, Vineberg K, Damen E, Van As CJ, et al. Dysphagia and aspiration after chemoradiotherapy for head-and-neck cancer: which anatomic structures are affected and can they be spared by IMRT? *Int J Radiat Oncol Biol Phys*. 2004;60(5):1425-39.

96. Hind JA, Nicosia MA, Roecker EB, Carnes ML, Robbins J. Comparison of effortful and noneffortful swallows in healthy middle-aged and older adults. *Arch Phys Med Rehabil*. 2001;82(12):1661-5.

97. Fujii M, Logemann J. Effect of a tongue-holding maneuver on posterior pharyngeal wall movement during deglutition. *Am J Speech Lang Pathol*. 1995;43:389-95.

98. Bandy WD, Irion JM, Briggler M. The effect of static stretch and dynamic range of motion training on the flexibility of the hamstring muscles. *J Orthop Sports Phys Ther*. 1998;27(4):295-300.

99. Bandy WD, Irion JM. The effect of time on static stretch on the flexibility of the hamstring muscles. *Phys Ther*. 1994;74(9):845-50.

100. Bandy WD, Irion JM, Briggler M. The effect of time and frequency of static stretching on flexibility of the hamstring muscles. *Phys Ther*. 1997;77(10):1090-6.

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101. Feland JB, Myrer JW, Schulthies SS, Fellingham GW, Measom GW. The effect of duration of stretching of the hamstring muscle group for increasing range of motion in people aged 65 years or older. *Phys Ther.* 2001;81(5):1110-7.

102. Michelotti A, Steenks MH, Farella M, Parisini F, Cimino R, Martina R. The additional value of a home physical therapy regimen versus patient education only for the treatment of myofascial pain of the jaw muscles: short-term results of a randomized clinical trial. *J Orofac Pain.* 2004;18(2):114-25.

103. Ucar M, Sarp U, Koca I, Eroglu S, Yetisgin A, Tutoglu A, et al. Effectiveness of a home exercise program in combination with ultrasound therapy for temporomandibular joint disorders. *J Phys Ther Sci.* 2014;26(12):1847-9.

104. Shaffer SM, Brismee JM, Sizer PS, Courtney CA. Temporomandibular disorders. Part 2: conservative management. *J Man Manip Ther.* 2014;22(1):13-23.

105. McNeely ML, Armijo OS, Magee DJ. A systematic review of the effectiveness of physical therapy interventions for temporomandibular disorders. *Phys Ther.* 2006;86(5):710- 25.

106. van Hinte G, Leijendekkers RA, Te Molder B, Jansen L, Bol C, Merckx MAW, et al. Reproducibility of measurements on physical performance in head and neck cancer survivors; measurements on maximum mouth opening, shoulder and neck function, upper and lower body strength, level of physical mobility, and walking ability. *PLoS One.* 2020;15(9):e0233271.

107. Kraaijenga SA, Hamming-Vrieze O, Verheijen S, Lamers E, van der Molen L, Hilgers FJ, et al. Radiation dose to the masseter and medial pterygoid muscle in relation to trismus after chemoradiotherapy for advanced head and neck cancer. *Head Neck.* 2021 Jul;29(7):3793-3800.

108. Sandler ML, Lazarus CL, Ru M, Sharif KF, Yue LE, Griffin MJ, et al. Effects of jaw exercise intervention timing on outcomes following oral and oropharyngeal cancer surgery: Pilot study. *Head Neck.* 2019 Nov;41(11):3806-3817.

109. Kumari S, Reddy DC, Paul ST. The normal range of maximal incisal opening in pediatric population and its association with physical variables. *Ann Afr Med.* 2019;18(3):153-7.

110. Cuellar ME, Oommen E. Objective physiological measures of lingual and jaw function in healthy individuals and individuals with dysphagia due to neurodegenerative diseases. *MethodsX.* 2021;8:101461.

111. McMillan H, Barbon CEA, Cardoso R, Sedory A, Buoy S, Porsche C, et al. Manual Therapy for Patients With Radiation-Associated Trismus After Head and Neck Cancer. *JAMA Otolaryngol Head Neck Surg.* 2022 May 1;148(5):418-425.

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<b>Autores</b>	<b>Título</b>	<b>Produto</b>	<b>Comentário</b>
Kumari et al., 2019 (109)	The normal range of maximal incisal opening in pediatric population and its association with physical variables	ROM scale	Observational prospective cohort study [Estudo de coorte prospectivo observacional]
Sandler et al., 2019 (108)	Effects of jaw exercise intervention timing on outcomes following oral and oropharyngeal cancer surgery: Pilot study	ROM scale	Prospective pilot study [Estudo piloto prospectivo]
Ezzat et al., 2021 (37)	The role of exercise therapy in managing postradiotherapy trismus in head and neck cancer	TheraBite® Jaw Mobilizer, ROM scale	Trismus rehabilitation, case report [Reabilitação do trismo: relato de caso]
Montalvo et al., 2020 (56)	Impact of exercise with TheraBite® device on trismus and health-related quality of life: A prospective study	TheraBite® Jaw Mobilizer	Trismus rehabilitation, prospective study [Reabilitação do trismo: estudo prospectivo]
Shao et al., 2020 (86)	Exercise therapy for cancer treatment- induced trismus in patients with head and neck cancer: A systematic review and meta-analysis of randomized controlled trials	TheraBite® Jaw Mobilizer	Trismus rehabilitation, systematic review/meta-analysis [Reabilitação do trismo: revisão sistemática/meta-análise]
Lloyd Morris et al., 2020 (74)	Improving Temporomandibular Range of Motion in People With Duchenne Muscular Dystrophy and Spinal Muscular Atrophy	TheraBite® Jaw Mobilizer, ROM scale	Rehabilitation of improved ROM, case series [Reabilitação com melhora da amplitude de movimento, série de casos]

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Reed et al., 2020 (65)	Dental management of scleroderma patients using pentoxifylline plus vitamin E with and without TheraBite® to reduce trismus: Two case reports and brief review of literature	TheraBite® Jaw Mobilizer	Trismus rehabilitation, case series [Reabilitação do trismo, série de casos]
Karsten et al., 2020 (89)	Long-term swallowing, trismus, and speech outcomes after combined chemoradiotherapy and preventive rehabilitation for head and neck cancer	TheraBite® Jaw Mobilizer, ROM scale	Trismus prevention and rehabilitation, prospective study [Prevenção e reabilitação do trismo: estudo prospectivo]
Van der Geer et al., 2020 (43)	The use of stretching devices for treatment of trismus in head and neck cancer patients: a randomized controlled trial	TheraBite® Jaw Mobilizer	Trismus rehabilitation, randomized controlled trial [Reabilitação do trismo, ensaio clínico randomizado controlado]
Bragante et al., 2020 (58)	Efficacy of exercise therapy during radiotherapy to prevent reduction in mouth opening in patients with head and neck cancer: A randomized controlled trial	TheraBite® Jaw Mobilizer	Trismus prevention, case study [Prevenção do trismo: estudo de caso]
Van Hinte et al., 2020 (106)	Reproducibility of measurements on physical performance in head and neck cancer survivors; measurements on maximum mouth opening, shoulder and neck function, upper and lower body strength, level of physical mobility, and walking ability	ROM scale	Reproducibility measurements on physical performance in HNC survivors [Medidas de reprodutibilidade do desempenho físico em sobreviventes de câncer de cabeça e pescoço]
Cuellar & Oommen, 2021 (110)	Objective physiological measures of lingual and jaw function in healthy individuals and individuals with dysphagia due to neurodegenerative diseases	ROM scale	Dysphagia treatment and jaw movement, clinical protocol [Tratamento da disfagia e movimento da mandíbula: protocolo clínico]
Warwas et al., 2021 (78)	Trismus Due to Hyperplastic Coronoid Process: Series of Five Cases and Review of Literature	TheraBite® Jaw Mobilizer	MIO rehabilitation, Case series [Reabilitação MIO, Série de casos]

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Chee et al., 2021 (47)	Interventions for Trismus in Head and Neck Cancer Patients: A Systematic Review of Randomized Controlled Trials	TheraBite® Jaw Mobilizer	Trismus rehabilitation, systematic review [Reabilitação do trismo: revisão sistemática]
Sydow et al., 2021 (75)	A mixed method study exploring patient satisfaction and feasibility of two exercise programmes in systemic sclerosis-associated microstomia	TheraBite® Jaw Mobilizer	Trismus rehabilitation, randomized clinical trial [Reabilitação do trismo: ensaio clínico randomizado]
Karlsson et al., 2021 (93)	Jaw exercise therapy for the treatment of trismus in head and neck Cancer: a prospective three-year follow-up study	TheraBite® Jaw Mobilizer	Trismus rehabilitation and prevention, prospective study [Reabilitação e prevenção do trismo: estudo prospectivo]
Aryanpour et al., 2022 (79)	Mandibular Myofibroma and Severe Trismus: A Complex Case and Review of Complications	TheraBite® Jaw Mobilizer	Pediatric trismus rehabilitation, case report [Reabilitação de trismo pediátrico: relato de caso]
Wang et al., 2022 (85)	Exercise for Trismus Prevention in Patients with Head and Neck Cancer: A Network Meta-Analysis of Randomized Controlled Trials	TheraBite® Jaw Mobilizer	Trismus prevention, meta-analysis [Prevenção do trismo, meta-análise]
McMillan et al., 2022 (111)	Manual Therapy for Patients With Radiation-Associated Trismus After Head and Neck Cancer	ROM scale	Trismus rehabilitation, retrospective case series [Reabilitação do trismo: série de casos retrospectiva]
Chambers et al., 2022 (94)	Trismus therapy devices: A systematic review	TheraBite® Jaw Mobilizer	Trismus prevention and rehabilitation, literature review [Prevenção e reabilitação do trismo: revisão da literatura]
Yim et al., 2022 (66)	Treatment of Temporomandibular Joint Heterotopic Ossification: A Novel Protocol With Multimodal Therapy Based on Literature Review and Presentation of a Unique Case Report	TheraBite® Jaw Mobilizer	TMJ heterotopic ossification protocol [Protocolo de ossificação heterotópica da ATM]
Karsten et al., 2022 (84)	Dysphagia, trismus and speech impairment following radiation-based treatment for advanced stage oropharyngeal carcinoma: a one-year prospective evaluation	TheraBite® Jaw Mobilizer	Trismus prevention, prospective study/Literature review [Prevenção do trismo: estudo prospectivo/Revisão da literatura]

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© 2023/06 Coloplast A/S. Todos os direitos reservados. Atos e o logotipo da Atos Medical são marcas registradas da Coloplast A/S.
Atos Medical AB, Hyllie Boulevard 1 7, SE-215 32 Malmö, Suécia
Atos Medical AB, Kraftgatan 8, SE-242 35 Hörby, Suécia Tel: +46 (0) 415 198 00